



Nitrile exam gloves with the quality you require, the performance and protection you'll love, and the value you expect.

Our Syntrile® XTS Powder Free 3.5g Nitrile Chemo Exam Gloves provide enhanced protection against chemotherapy drugs* -- yet are thin for greater comfort and conformability. They feature a textured finish at the fingertip to promote slip resistance in wet or dry conditions.

Cypress Medical Products LLC is committed to offering exceptional value to customers by delivering high-quality, competitively-priced products, which have been put through rigorous quality testing. Our factory audits meet the FDA's cGMP guidelines. Our gloves meet or exceed ASTM standards.

Features

- Powder free
- Beaded cuff
- Super soft
- Textured finish
- Single use/disposable

Benefits

- Alleviates powder irritation concerns
- Extra strength and ease in glove donning
- Enhanced comfort
- Provides a better grip, wet or dry, and improves tactile sensitivity and dexterity
- Eliminates cross-contamination risk
- Minimizes latex allergy and protein concerns



Syntrile™ XTS Powder Free Nitrile Chemo Exam Gloves Specifications

Fingertip Thickness	.07mm (2.76 mil)
Palm Thickness	.05mm (1.97 mil)
Cuff Thickness	.05mm (1.97 mil)
Color	Azure Blue
Textured	Fingertips

10
boxes per case

Small.....	27-12	200/box
Medium ...	27-14	200/box
Large.....	27-16	200/box
X-Large.....	27-18	180/box

For more information, please contact your Cypress Customer Service Representative at 800.334.3646 or your local distributor.

CAUTION: THESE GLOVES MAY PUNCTURE OR TEAR IF DONNED IMPROPERLY OR BY BEING CONTACT WITH SHARP FINGERNAILS OR OBJECTS SUCH AS FINGER RINGS, ETC.

CAUTION: ALTHOUGH EXAMINATION GLOVES PROVIDE AN EXCELLENT PROTECTIVE BARRIER AGAINST BLOOD BORNE PATHOGENS. THEY ARE NOT INTENDED FOR PROLONGED, DIRECT EXPOSURE TO HARSH CHEMICALS, WHERE HEAVY DUTY OR INDUSTRIAL GLOVES ARE REQUIRED.

CAUTION: COMPONENTS USED IN GLOVE MANUFACTURING MAY CAUSE ALLERGIC REACTIONS IN SOME USERS.

* Tested and approved for use as chemo exam gloves. See Reverse side for chemotherapy test data

Syntrile™ XTS Nitrile Powder Free Chemo Exam Glove Specifications

Dimensions <i>(millimeters)</i>	Length	SKU #	240mm (min.)
	Small Width	27-12	87mm
	Medium Width	27-14	97mm
	Large Width	27-16	107mm
	X-Large Width	27-18	115mm
	Fingertip Thickness		.07mm (2.76 mil)
	Palm Thickness		.05mm (1.97 mil)
	Cuff Thickness		.05mm (1.97 mil)
Strength	Pinhole AQL (water leak test)		2.5
	Tensile Strength Min(before aging)		18 Mpa
	Tensile Strength Min(after aging)		16 Mpa
Elasticity	Elongation Min(before aging)		500%
	Elongation Min(after aging)		400%
Color			Azure Blue
Textured			Fingertips

Chemotherapy Drug Permeation Resistance

(average normalized breakthrough time in minutes)

Carmustine (BCNU) (3,300ppm or 3.3mg/mL)	NR*
Cisplatin (1,000ppm or 1.0mg/mL)	>240
Cyclophosphamide (Cytoxan) (20,000ppm or 20.0mg/mL)	>240
Dacarbazine (DTIC) (10,000ppm or 10.0mg/mL)	>240
Doxorubicin Hydr ochloride (2,000ppm or 2.0mg/mL)	>240
Etoposide (Toposar) (20,000ppm or 20.0mg/mL)	>240
5-Fluor ouracil (50,000ppm or 50.0mg/mL)	>240
Paclitaxel (Taxol) (6,000ppm or 6.0mg/mL)	>240
Thiotepa (10,000ppm or 10.0mg/mL)	NR*

Tested for use with Chemotherapy Drugs using ASTM D6978-05. Gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of drugs being used.

Chemical Permeation Data

Please review the drug labeling or Material Safety Data Sheet for the chemotherapy drugs to be used to determine the required level of protection for the intended use.

Meets or Exceeds ASTM D-6319

***Not recommended for use with Thiotepa and Carmustine.**

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by Cypress Medical Products LLC or its affiliates ("Cypress"). Cypress is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from Cypress's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and Cypress disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by Cypress. All trademarks and registered trademarks are the property of their respective owners.

Cypress Medical Products LLC
1400 South Wolf Road
Suite 200
Wheeling, IL 60090

800-334-3646
sales@cypressmed.com
cypressmed.com
© 2013 Cypress Medical Products LLC